




Whatever you want, You will find it!

Sorenson™ Filter Tips



Standard Filter Tips

Low Binding Filter Tips

OneTouch Filter Tips

Sorenson Filter Tips

Sorenson사는 높은 품질의 Pipette tip을 공급하고 있으며, 특히 Filter Tip 산업에서 탁월하고 혁신적인 제품을 공급하고 있습니다.
다양한 Filter Tip을 실험목적에 맞추어, 합리적인 가격으로 경험해 보십시오.

필터팁 적용 분야: 일반적인 바이오 실험, PCR, 바이러스 검출, 질병 진단, 방사선 동위원소 실험, 오염방지, 시료확보, 피펫 오염 방지 등 다양한 분야

Sorenson Filter Tip의 종류

Filter Tip 종류	주요 특징
Standard Filter Tips	일반적인 실험에 사용(합리적인 가격)
Low Binding Filter Tips	Sorenson사의 주요 특징인 Low Binding Technology를 Filter Tip에 적용시켜 DNA, Peptide, Protein, Cell 관련 실험 시 시료가 Tip에 묻어나는 정도를 크게 낮추어 줍니다.
OneTouch Filter Tips	피펫과 결합되는 부분을 TPE로 만든 팁으로서 가벼운 힘으로도 완벽한 결합이 됩니다. 8-ch, 12-ch 피펫 사용시 등 실험량이 많은 경우, 정확하고 정확한 볼륨을 원하실 경우 사용됩니다.

Sorenson Filter Tip의 특징

▷ 다양한 필터팁 제형, 제품군 및 높은 피펫 호환성

▶ 필터가 팁 상단에 위치하여 오버피펫팅 방지(250ul, 1250ul 등)

▷ Locking Ring: 필터를 단단하게 고정함

▶ Cellulose 미첨가 필터(Cellulose : PCR Inhibitor로 작용)

** 학술자료 : <http://www.sorbio.com/performance.cfm>

(위 url을 참조하시면 탄저균, 방사선 오염에 대한 효과적인 필터효과에 대한 자료를 보실 수 있습니다.)

▷ CellRack Package

▶ 타사보다 얇은 필터 두께: 필터에 샘플이 닿는 것을 방지

Sorenson Certified Free of



Sorenson 제품은 최고의 Certification을 제공합니다. 최고 품질의 제품으로 연구에 안심을 더하십시오. 다음의 5가지 항목에 대해 Free 합니다.

- ▶ RNase Free
- ▶ DNase Free
- ▶ Human DNA Free
- ▶ PCR Inhibitors Free
- ▶ Pyrogen Free

Sorenson Filter Tips - Standard

Cat. No.	Description	Vol. Range	Graduations	Rack Type
38000	10ul LongReach Filter Tip	0.1-10ul	2/5/10ul	Normal
30340	10ul Precision Filter Tip	0.1-10ul	1/2/5/10ul	Hinged
14210	20ul Filter Tip	1-20ul	No Graduations	Normal
36060	100ul Filter Tip	1-100ul	No Graduations	Normal
30550	200ul NX Filter Tip	1-250ul	10/50/100ul	Normal
17370	Multiguard 3 200ul Tip (20, 50, 100, 200ul Pipettor에 모두 사용 가능)	1-200ul	10/50/100/200ul	Normal
30510-X	300ul XT Filter Tip	1-300ul	10/50/100/200ul	Normal
14200	1000ul Filter Tip	100-1000ul	No Graduations	Normal
34000	1000ul XT Filter Tip	50-1250ul	200/500/1000ul	Hinged

Sorenson Filter Tips - Low Binding










Cat. No.	Description	Vol. Range	Graduations	Rack Type
38000T	10ul LongReach Filter Tip	0.1-10ul	2/5/10ul	Normal
14210T	20ul Filter Tip	1-20ul	No Graduations	Normal
36060T	100ul Filter Tip	1-100ul	No Graduations	Normal
30550T	200ul NX Filter Tip	1-250ul	10/50/100ul	Normal
17370T	Multiguard 3 200ul Tip (20, 50, 100, 200ul Pipettor에 모두 사용 가능)	1-200ul	10/50/100/200ul	Normal
30510T-X	300ul XT Filter Tip	1-300ul	10/50/100/200ul	Normal
14200T	1000ul Filter Tip	100-1000ul	No Graduations	Normal

Sorenson Filter Tips - OneTouch





Cat. No.	Description	Vol. Range	Graduations	Rack Type
78510	3ul OneTouch Filter Tip	0.1-3ul	1ul	Normal
10320	10ul OneTouch Filter Tip	1-20ul	1/5/10ul	Normal
10350	200ul OneTouch Filter Tip	1-330ul	10/50/100ul	Normal
10380	1000ul OneTouch Filter Tip	50-1250ul	200/500/1000ul	Hinged

Pipetting like you've never felt it

Sorenson Filter Tips - Standard, Low Binding

	10ul LongReach	38000 38000T	45mm
	10ul Precision	30340	45mm
	20ul	14210 14210T	51mm
	100ul	36060 36060T	51mm
	200ul	30550 30550T	59mm
	Multiguard 3 200ul	17370 17370T	78mm
	300ul	30510-X 30510T-X	78mm
	1000ul	14200 14200T	72mm
	1000ul XT	34000	102mm

Sorenson Filter Tips - OneTouch

	3ul	78510	31mm
	10ul	10320	45mm
	200ul	10350	59mm
	1000ul	10380	102mm

[서울사무실] (04799) 서울시 성동구 광나루로6길 35 우림이비즈센터 209호

[양주물류센터] (11423) 경기도 양주시 광적면 부흥로 521-57

Tel. 02-6409-9338 / Fax. 02-6499-9338 / E-mail. info@b2bio.co.kr

www.b2bio.co.kr / blog.naver.com/b2bioinc